# PMH



# FLEXIVAC® PRODUCTS CATALOG



MH-03P-02F REV.01

# Care with Excellence

PMH – Produtos Médicos Hospitalares, S.A. | +351 263 650 680 | www.pmh.pt | info@pmh.pt



#### SUCTION BAGS

**FLEXIVAC® SYSTEMS** 

# ENSURE ASEPTIC AND SAFE HANDLING DURING SUCTION OF ORGANIC FLUIDS.

**FLEXIVAC** IS A CLOSED, FLEXIBLE, SAFE SUCTION SYSTEM FOR ASPIRATION OF ORGANIC FLUIDS AND FUMES. PREVENTS CROSS-CONTAMINATION.

#### **GREATER SECURITY**

- IT IS A CLOSED SYSTEM.
- SECURE CONNECTION SYSTEM.
- PREVENTION OF CROSS-CONTAMINATION.
- CHECK VALVE AT PATIENT CONNECTION.
- ANTIBACTERIAL FILTER.
- HYDROPHOBIC FILTER.

#### USE

- EASY INSTALLATION AND PLACEMENT OF BAGS.
- HANDLES FOR TRANSPORTATION AND REMOVAL.
- THE BELLOWS SYSTEM ALLOWS AUTO-DISCONNECTION WHEN THE BAG IS FULL.
- CLOSED SYSTEM, EASY REMOVAL WITHOUT THE NEED TO PUT ON CAPS, AVOIDING ACCIDENTAL SPILLS.

#### VARIOUS CAPABILITIES

- WITH FLEXIBLE BAGS FROM 1L TO 3L.
- VERSIONS WITH HOLES FOR INTRODUCING SOLIDIFYING POWDER AND/OR COLLECTING SAMPLES.
- INCORPORATION OF ADAPTER FOR DIFFERENT DIAMETERS OF THE SUCTION TUBE, TO REDUCE THE POSSIBILITY OF BENDING AND FOR BETTER POSITIONING.





# SUCTION BAGS

#### 704.2000 - FLEXIBLE BAG FOR ASPIRATION OF ORGANIC FLUIDS WITH LIQUID STOP FILTER

- FLEXIBLE BAG FOR COLLECTION OF FLUIDS
- GREEN FLEXIVAC CURVE
- FLEXIVAC CAP WITHOUT HOLE, GREEN
- LIQUID STOP FILTER
- ANTI-FOAM
- UNITS PER BOX: 70

# 704.2001 - FLEXIBLE BAG FOR ASPIRATION OF ORGANIC FLUIDS WITH HOLE AND LIQUID STOP FILTER

- FLEXIBLE BAG FOR COLLECTION OF FLUIDS
- GREEN FLEXIVAC CURVE
- FLEXIVAC CAP WITH HOLE, GREEN
- LIQUID STOP FILTER
- ANTI-FOAM
- UNITS PER BOX: 70

# 704.1861 - BAG FOR REABSORPTION OF TISSUE AND FAT

- FLEXIBLE BAG STERILE SUITABLE FOR REABORTION OF
- TISSUE AND FAT FROM PLASTIC SURGERIES
- GREEN FLEXIVAC CURVE
- FLEXIVAC CAP WITH HOLE, GREEN
- LIQUID STOP FILTER
- UNITS PER BOX: 60









# SUCTION BAGS

#### **3L CAPACITY**

## 704.3000 - FLEXIBLE BAG FOR ASPIRATION OF ORGANIC FLUIDS WITH LIQUID STOP FILTER

- FLEXIBLE BAG FOR COLLECTION OF FLUIDS
- GREEN FLEXIVAC CURVE
- FLEXIVAC CAP WITHOUT HOLE, GREEN
- LIQUID STOP FILTER
- ANTI-FOAM
- UNITS PER BOX: 50

# 704.3001 - FLEXIBLE BAG FOR ASPIRATION OF ORGANIC FLUIDS WITH HOLE AND LIQUID STOP FILTER

- FLEXIBLE BAG FOR COLLECTION OF FLUIDS
- GREEN FLEXIVAC CURVE
- FLEXIVAC CAP WITH HOLE, GREEN
- LIQUID STOP FILTER
- ANTI-FOAM
- UNITS PER BOX: 50







#### **TANDEM SYSTEM**


PRODUCTS WITH SERIES CONNECTION SYSTEM

ENSURE ASEPTIC AND SAFE HANDLING DURING SUCTION OF ORGANIC FLUIDS.

PMH SUCTION BAGS, WITH SERIES CONNECTION, ARE USED TO SUCTION ORGANIC FLUIDS, IN A SAFE AND EFFICIENT WAY, ALLOWING ASPIRATION.

BEING CONNECTED IN SERIES, FROM BAG TO BAG, THERE IS NO NEED FOR DISCONNECTION.







#### SUCTION BAGS 1L / 1,5L / 2L CAPACITY

## 704.2010 - FLEXIBLE BAG FOR ASPIRATION OF ORGANIC FLUIDS, FOR SERIES USE ONLY

- FLEXIBLE BAG FOR COLLECTION OF FLUIDS
- FLEXIVAC CAP WITHOUT HOLE, WHITE
- ANTI-FOAM

• 23CM TRANSPARENT SUCTION TUBE, WITH CLAMP AND GREEN FLEXIVAC CURVED CONNECTORS • UNITS PER BOX: 50

# 704.2011 - FLEXIBLE BAG FOR ASPIRATION OF ORGANIC FLUIDS WITH HOLE, FOR SERIES USE ONLY

- FLEXIBLE BAG FOR COLLECTION OF FLUIDS
- FLEXIVAC CAP WITH HOLE, WHITE
- ANTI-FOAM
- 23CM TRANSPARENT SUCTION TUBE, WITH CLAMP AND GREEN FLEXIVAC CURVED CONNECTORS
- UNITS PER BOX: 50





## SUCTION BAGS

#### **3L CAPACITY**

## 704.3010 - FLEXIBLE BAG FOR ASPIRATION OF ORGANIC FLUIDS, FOR SERIES USE ONLY

- FLEXIBLE BAG FOR COLLECTION OF FLUIDS
- FLEXIVAC CAP WITHOUT HOLE, WHITE
- ANTI-FOAM
- 23CM TRANSPARENT SUCTION TUBE, WITH CLAMP AND GREEN FLEXIVAC CURVED CONNECTORS
- UNITS PER BOX: 50





**FLEXIVAC® SYSTEMS** 

# ENSURE ASEPTIC AND SAFE HANDLING DURING SUCTION OF ORGANIC FLUIDS.

FLEXIVAC CONTAINERS ARE REUSABLE, HIGHLY RESISTANT AND EQUIPPED WITH A VISIBLE SCALE. AVAILABLE IN CAPACITIES OF 1L, 1.5L, 2L AND 3L, THEY ARE IDEAL FOR VARIOUS VACUUMING NEEDS, INCLUDING IN CLOSED ENVIRONMENTS SUCH AS AMBULANCES AND INTENSIVE CARE ROOMS.

THEY OFFER THE OPTION WITH A TAP TO CLOSE THE VACUUM, MAKING INSTALLATION AND REMOVAL EASIER. FLEXIVAC SYSTEMS PROVIDE EASE, COMFORT AND SAFETY, AVOIDING CHANGING BAGS DURING SURGICAL PROCEDURES AND ADAPTING TO THE INDIVIDUAL REQUIREMENTS OF EACH SUCTION SITUATION.

THEY ARE USED TO EQUIP FLEXIVAC ROTATED SUPPORTS, WHICH HOLD UP TO 8 CUPS OF 3L EACH, ADAPTED TO THE SPECIFIC NEEDS OF EACH USER. WITH A SIMPLE AND ROBUST DESIGN, THEY GUARANTEE SAFE ASPIRATION OF SECRETIONS, MAKING PMH SUCTION BAGS EASIER TO HANDLE.

VARIOUS CONFIGURATIONS AVAILABLE TO MEET THE SPECIFIC NEEDS OF EACH PROCEDURE.







#### **1L CONTAINERS**

#### 705.0874 - FLEXIVAC® 1L CONTAINER

• 1L FLEXIVAC CONTAINER

• MOUNTING BRACKET



- 1L FLEXIVAC CONTAINER
- MOUNTING BRACKET
- TAP

## 705.1495 - FLEXIVAC® 1L CONTAINER WITH TUBE

- 1L FLEXIVAC CONTAINER
- MOUNTING BRACKET
- 34CM GREEN VACUUM TUBE, WITH CONNECTOR

# 705.1650- FLEXIVAC® 1L CONTAINER WITH TUBE AND TAP

- 1L FLEXIVAC CONTAINER
- MOUNTING BRACKET
- TAP
- 34CM GREEN VACUUM TUBE, WITH CONNECTOR











#### 1,5L CONTAINERS

#### 705.0875 - FLEXIVAC® 1,5L CONTAINER

- 1,5L FLEXIVAC CONTAINER
- MOUNTING BRACKET

## 705.2284 - FLEXIVAC® 1,5L CONTAINER WITH TAP


- 1,5L FLEXIVAC CONTAINER
- MOUNTING BRACKET
- TAP


# 705.1496 - FLEXIVAC® 1,5L CONTAINER WITH TUBE


- 1,5L FLEXIVAC CONTAINER
- MOUNTING BRACKET
- 34CM GREEN VACUUM TUBE, WITH CONNECTOR


# 705.1651 - FLEXIVAC® 1,5L CONTAINER WITH TUBE AND TAP

- 1,5L FLEXIVAC CONTAINER
- MOUNTING BRACKET
- TAP
- 34CM GREEN VACUUM TUBE, WITH CONNECTOR











#### **2L CONTAINERS**

#### 705.1762 - FLEXIVAC® 2L CONTAINER

• 2L FLEXIVAC CONTAINER

MOUNTING BRACKET

# 705.2285 - FLEXIVAC® 2L CONTAINER WITH TAP

- 2L FLEXIVAC CONTAINER
- MOUNTING BRACKET
- TAP

# 705.1763 - FLEXIVAC® 2L CONTAINER WITH TUBE

- 2L FLEXIVAC CONTAINER
- MOUNTING BRACKET
- 34CM GREEN VACUUM TUBE, WITH CONNECTOR

# 705.1592 - FLEXIVAC® 2L CONTAINER WITH TUBE AND TAP

- 2L FLEXIVAC CONTAINER
- MOUNTING BRACKET
- TAP
- 34CM GREEN VACUUM TUBE, WITH CONNECTOR











#### **3L CONTAINERS**

705.0876 - FLEXIVAC® 3L CONTAINER

• 3L FLEXIVAC CONTAINER

• MOUNTING BRACKET

## 705.2286 - FLEXIVAC® 3L CONTAINER WITH TAP

• 3L FLEXIVAC CONTAINER

- MOUNTING BRACKET
- TAP

## 705.1392 - FLEXIVAC® 3L CONTAINER WITH TUBE

- 3L FLEXIVAC CONTAINER
- MOUNTING BRACKET
- 34CM GREEN VACUUM TUBE, WITH CONNECTOR

# 705.1652 - FLEXIVAC® 3L CONTAINER WITH TUBE AND TAP

- 3L FLEXIVAC CONTAINER
- MOUNTING BRACKET
- TAP
- 34CM GREEN VACUUM TUBE, WITH CONNECTOR





STROLLER SUPPORTS FLEXIVAC®

FLEXIVAC® SYSTEMS

705.0877 - ROTATED SUPPORT WITH CAPACITY FOR 2 FLEXIVAC® CONTAINERS



705.0878- ROTATED SUPPORT WITH CAPACITY FOR 4 FLEXIVAC® CONTAINERS

705.0879 - ROTATED SUPPORT WITH CAPACITY FOR 6 FLEXIVAC® CONTAINERS

705.0880 - ROTATED SUPPORT WITH CAPACITY FOR 8 FLEXIVAC® CONTAINERS







12



#### **FLEXIVAC® ACCESSORIES**

#### FLEXIVAC® SYSTEMS

#### 705.1872 - EXTERNAL FILTER TO PROTECT THE SUCTION LINE - LINEAVAC

- FLEXIVAC® FILTER PROTECTION BODY AND TAP
- LARGE LIQUID STOP FILTER
- UNITS PER BOX: 50



#### 705.4176 - FRAGMENT FILTER - SOLIDFLEX

- FLEXIVAC® FILTER PROTECTION BODY AND TAP
- HARVEST FILTER
- SUCTION TUBE CONNECTOR
- UNITS PER BOX: 50



#### 705.4478 - SOLIDIFYING POWDER IN 20G BAG

• WATER-SOLUBLE BAG WITH 20G OF SOLIDIFYING POWDER • UNITS PER BOX: 100

#### 705.2015 - SOLIDIFYING POWDER IN 30G BAG

• WATER-SOLUBLE BAG WITH 30G OF SOLIDIFYING POWDER

• UNITS PER BOX: 100







#### **FLEXIVAC® ACCESSORIES**

#### FLEXIVAC® SYSTEMS

#### 705.0884 - FLEXIVAC® WALL SUPPORT

• UNITS PER BOX: ACCORDING TO CUSTOMER REQUEST



• UNITS PER BOX: ACCORDING TO CUSTOMER REQUEST



