

NEOSEGUR® PRODUCTS

PRODUCTS CATALOG

Care with Excellence



NEOSEGUR® EXTENSIONS

NEOSEGUR® SYSTEMS

**ENSURE SAFE MEDICATION ADMINISTRATION,
WITHOUT MIXING MEDICATIONS IN THE LINES.**

- THE LINES COME TOGETHER AT THE LUER LOCK THAT CONNECTS THE CATHETER
- THE ADMINISTERED DRUGS DO NOT MIX DURING THE PROLONGATION
- DOES NOT ROTATE, MAKES IT EASY TO PLACE WITH JUST ONE HAND
- ALLOWS THE MAINTENANCE OF THE SELECTED FLOW RATE WITHOUT RETURN
- COLOR-CODED LINES ON THE LUER LOCK
- OPTIMIZED IDENTIFICATION OF DRUGS TO BE ADMINISTERED IN THE LINES AND POSSIBILITY OF REMOVAL FOR GREATER CONVENIENCE
- MOBILE LUER LOCK WITH PURGE FILTER
- COLOR-CODED OPEN CLAMPS
- SAFETY AND COMFORT



NEOSEGUR®

1-WAY

**306.1111 - 1-WAY NEOSEGUR® EXTENSION,
WITH BIDIRECTIONAL VALVE, MOVABLE
LUER LOCK AND PURGE FILTER**

- PROTECTIVE CAP WITH PURGE FILTER
- DOUBLE POE/PVC TUBE 1.00X3.00MM, 15CM
- WHITE CLAMP
- NEUTROVAL® BIDIRECTIONAL VALVE
- UNITS PER BOX: 150





NEOSEGUR®

2-WAY

306.2828 - 2-WAY NEOSEGUR® EXTENSION, WITH VAR AND MOBILE LUER LOCK AND PURGE FILTER

- PROTECTIVE CAP WITH PURGE FILTER
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- WHITE AND BLUE CLAMP
- NON-RETURN VALVE (VAR)
- VALVE PROTECTION CAP
- UNITS PER BOX: 150




306.2862 - 2-WAY NEOSEGUR® EXTENSION, WITH BIDIRECTIONAL VALVE AND MOVABLE LUER LOCK AND PURGE FILTER

- PROTECTIVE CAP WITH PURGE FILTER
- POE/PVC DOUBLE TUBE 1.50X3.00MM, 10CM
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- RED AND BLUE CLAMP
- NEUTROVAL® BIDIRECTIONAL VALVE
- UNITS PER BOX: 150



306.2820 - 2-WAY NEOSEGUR® EXTENSION, ONE HIGH-FLOW EXTENSION, WITH BIDIRECTIONAL VALVE AND MOVABLE LUER LOCK AND PURGE FILTER

- PROTECTIVE CAP WITH PURGE FILTER
- THIN TUBE WITH BLUE LINE 1.80X3.00MM, 10CM
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- RED AND BLUE CLAMP
- NEUTROVAL® BIDIRECTIONAL VALVE
- UNITS PER BOX: 150



306.4465 - 2-WAY ORANGE OPAQUE NEOSEGUR® EXTENSION, WITH VAR AND BIDIRECTIONAL VALVE, MOVABLE LUER LOCK AND PURGE FILTER

- PROTECTIVE CAP WITH PURGE FILTER
- ORANGE POE/PVC DOUBLE TUBE 1.50X3.00MM, 10CM
- ORANGE POE/PVC DOUBLE TUBE 0.50X2.00MM, 10CM
- RED AND BLUE CLAMP
- MICROCLAVE® BIDIRECTIONAL VALVE
- NON-RETURN VALVE (VAR)
- VALVE PROTECTION CAP
- UNITS PER BOX: 150






NEOSEGUR®

3-WAY


306.2787 - 3-WAY NEOSEGUR® EXTENSION, BLUE LINE, WITH BIDIRECTIONAL VALVE, MOVABLE LUER LOCK AND PURGE FILTER

- PROTECTIVE CAP WITH PURGE FILTER
- THIN TUBE WITH BLUE LINE 1.80X3.00MM, 10CM
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- WHITE, BLUE AND RED CLAMP
- NEUTROVAL® BIDIRECTIONAL VALVE
- UNITS PER BOX: 150




306.2784 - 3-WAY NEOSEGUR® EXTENSION WITH 3 VAR, MOVABLE LUER LOCK AND PURGE FILTER

- PROTECTIVE CAP WITH PURGE FILTER
- THIN TUBE WITH BLUE LINE 1.80X3.00MM, 10CM
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- WHITE, BLUE AND RED CLAMP
- NON-RETURN VALVE (VAR)
- VALVE PROTECTION CAP
- UNITS PER BOX: 150




306.3928 - 3-WAY NEOSEGUR® EXTENSION WITH 2 VAR, MOVABLE LUER LOCK AND PURGE FILTER

- PROTECTIVE CAP WITH PURGE FILTER
- THIN TUBE WITH BLUE LINE 1.80X3.00MM, 10CM
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- BLUE CLAMP
- NON-RETURN VALVE (VAR)
- VALVE PROTECTION CAP
- UNITS PER BOX: 150



306.2786 - 3-WAY NEOSEGUR® EXTENSION, WITH 2 BIDIRECTIONAL VALVES AND 1 VAR, WITH MOVABLE LUER LOCK AND PURGE FILTER

- PROTECTIVE CAP WITH PURGE FILTER
- THIN TUBE WITH BLUE LINE 1.80X3.00MM, 10CM
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- WHITE, BLUE AND RED CLAMP
- NEUTROVAL® BIDIRECTIONAL VALVE
- NON-RETURN VALVE (VAR)
- VALVE PROTECTION CAP
- UNITS PER BOX: 150






NEOSEGUR®

4-WAY


306.4068 - 4-WAY NEOSEGUR® EXTENSION, WITH 4 VAR AND MOBILE LUER LOCK WITH PURGE FILTER

- PROTECTIVE CAP WITH PURGE FILTER
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- WHITE, BLUE, RED AND GREEN CLAMP
- NON-RETURN VALVE (VAR)
- VALVE PROTECTION CAP
- UNITS PER BOX: 150




306.4135 - 4-WAY NEOSEGUR® EXTENSION, WITH BIDIRECTIONAL VALVE, MOVABLE LUER LOCK AND PURGE FILTER

- PROTECTIVE CAP WITH PURGE FILTER
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- WHITE, BLUE, RED AND GREEN CLAMP
- NEUTROVAL® BIDIRECTIONAL VALVE
- UNITS PER BOX: 150



306.4466 - 4-WAY NEOSEGUR® OPAQUE ORANGE EXTENSION, WITH 3 VAR AND A BIDIRECTIONAL VALVE, WITH MOVABLE LUER LOCK AND PURGE FILTER

- PROTECTIVE CAP WITH PURGE FILTER
- ORANGE POE/PVC DOUBLE TUBE 0.50X2.00MM, 150CM
- WHITE, BLUE, RED AND GREEN CLAMP
- MICROCLAVE® BIDIRECTIONAL VALVE
- NON-RETURN VALVE (VAR)
- VALVE PROTECTION CAP
- UNITS PER BOX: 100





NEOSEGUR®

5-WAY


306.4148 - 5-WAY NEOSEGUR® EXTENSION, WITH 5 VAR, WITH MOVABLE LUER LOCK AND PURGE FILTER

- PROTECTIVE CAP WITH PURGE FILTER
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- WHITE, BLUE, RED AND GREEN CLAMP
- NON-RETURN VALVE (VAR)
- VALVE PROTECTION CAP
- UNITS PER BOX: 150



306.4149 - 5-WAY NEOSEGUR® EXTENSION, WITH BIDIRECTIONAL VALVE, MOVABLE LUER LOCK AND PURGE FILTER

- PROTECTIVE CAP WITH PURGE FILTER
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- WHITE, BLUE, RED AND GREEN CLAMP
- NEUTROVAL® BIDIRECTIONAL VALVE
- UNITS PER BOX: 150






NEOSEGUR®

2-WAY - CARESAFE® VALVE


306.4499 - 2-WAY NEOSEGUR® EXTENSION, WITH CARESAFE®

- GREEN PROTECTIVE CAP WITH PURGE FILTER
- DOUBLE POE/PVC TUBE 1.50X3.00MM, 15CM
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 15CM
- BLUE AND WHITE CLAMP
- 2 CARESAFE® BIDIRECTIONAL VALVES
- UNITS PER BOX: 150



306.4732 - 2-WAY NEOSEGUR® EXTENSION, WITH CLAMP AND CARESAFE®

- GREEN PROTECTIVE CAP WITH PURGE FILTER
- POE/PVC DOUBLE TUBE 2.50X4.10MM, 10CM
- POE/PVC DOUBLE TUBE 1.50X3.00MM, 10CM
- BLUE AND WHITE CLAMP
- 2 CARESAFE® BIDIRECTIONAL VALVES
- UNITS PER BOX: 150




NEOSEGUR®

3-WAY - CARESAFE® VALVE

306.4539 - 3-WAY NEOSEGUR® EXTENSION, WITH CARESAFE®

- GREEN PROTECTIVE CAP WITH PURGE FILTER
- POE/PVC DOUBLE TUBE 1.80X3.00MM, 10CM
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- BLUE, WHITE AND RED CLAMP
- 3 CARESAFE® BIDIRECTIONAL VALVES
- UNITS PER BOX: 150






NEOSEGUR®

4-WAY - CARESAFE® VALVE

306.4557 - 4-WAY NEOSEGUR® EXTENSION, WITH CARESAFE®

- GREEN PROTECTIVE CAP WITH PURGE FILTER
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- BLUE, WHITE, RED AND GREEN CLAMP
- 4 CARESAFE® BIDIRECTIONAL VALVES
- UNITS PER BOX: 150




NEOSEGUR®

5-WAY - CARESAFE® VALVE

306.4545 - 5-WAY NEOSEGUR® EXTENSION, WITH CARESAFE®

- GREEN PROTECTIVE CAP WITH PURGE FILTER
- POE/PVC DOUBLE TUBE 0.90X2.00MM, 10CM
- BLUE, WHITE, RED AND GREEN CLAMP
- 5 CARESAFE® BIDIRECTIONAL VALVES
- UNITS PER BOX: 150





VALVES CHARACTERISTICS

NEOSEGUR® SYSTEMS

806.2897 - NEUTROVAL® VALVE

- SPLIT-SEPTUM TECHNOLOGY VALVE
- BEFORE USE, THIS MEMBRANE MUST BE DECONTAMINATED WITH 70% ALCOHOL OR CHLORHEXIDINE.
- INJECTION/ASPIRATION WITH POINT OF INFECTION, FOR USE WITHOUT A NEEDLE;
- BARRIER AGAINST MICROORGANISMS IN PROLONGED USE;
- SIMPLE AND EASILY ACCESSIBLE - LUER LOCK AND LUER SLIP;
- COMPATIBLE WITH ALCOHOL, LIPIDS, CHLORHEXIDINE AND CYTOSTATICS;
- COMPATIBLE WITH MRI AND BLOOD.
- PRESSURE TYPE PRESSURE: NEUTRAL
- FLOW RATE: 165ML/MIN
- RESIDUAL VALUE: <0.04ML
- MAXIMUM PRESSURE: 27.58 BAR
- NUMBER OF ACTIVATIONS: 7 DAYS / 700 ACTIVATIONS
- UNITS PER BOX: 200



806.2300 - MICROCLAVE® VALVE

- MICROCLAVE VALVE WITH SLIP-SEPTUM TECHNOLOGY
- BEFORE USE, THIS MEMBRANE MUST BE DECONTAMINATED WITH 70% ALCOHOL OR CHLORHEXIDINE.
- INJECTION/ASPIRATION WITH DISINFECTION POINT, FOR USE WITHOUT A NEEDLE;
- BARRIER AGAINST MICROORGANISMS IN PROLONGED USE;
- SIMPLE HANDLING AND EASY ACCESS - LUER LOCK AND LUER SLIP;
- COMPATIBLE WITH ALCOHOL, LIPIDS, CHLORHEXIDINE AND CYTOSTATICS;
- COMPATIBLE WITH MRI AND BLOOD.
- PRESSURE TYPE PRESSURE: NEUTRAL
- FLOW RATE: 165 ML/MIN
- RESIDUAL VOLUME: <0.04 ML
- MAX PRESSURE: 27.58 BAR
- NUMBER OF ACTIVATIONS: 7 DAYS / 700 ACTIVATIONS
- UNITS PER BOX: 200





VALVES CHARACTERISTICS

NEOSEGUR® SYSTEMS

806.0677 - ANTI-RETURN VALVE (VAR)

- NON-RETURN VALVE.
- LIPID RESISTANT;
- COMPATIBLE WITH CHEMOTHERAPY AND NMR.
- FLOW ≥ 150 ML/MIN.
- ≥ 90 ML/MIN. (GLUCOSE 40%)
- MAXIMUM PRESSURE RESISTANCE: 6BAR
- OPENING PRESSURE: ≤ 12 MBAR
- UNITS PER BOX: 200



806.4409 - CARESAFE® VALVE

- SPLIT-SEPTUM TECHNOLOGY VALVE. LUER ACTIVATION OR LUER LOCK WITH DIVIDED SEPTUM;
- CLOSED FLAT SEAL THAT PROVIDES A MICROBIAL BARRIER;
- THE VALVE IS OPEN DURING ASPIRATION OR DURING INFUSION, AND CLOSES WHEN IT IS SEPARATED FROM THE INFUSION DEVICE;
- SIMPLE HANDLING AND EASY ACCESS - LUER LOCK AND LUER SLIP;
- SUITABLE FOR TPN, BLOOD AND ONCOLOGY INFUSIONS;
- COMPATIBLE WITH LIPIDS, ALCOHOL, CHLORHEXIDINE AND CYTOSTATICS;
- MRI COMPATIBLE.
- PRESSURE TYPE PRESSURE: NEUTRAL
- FLOW RATE: 110-150ML/MIN
- RESIDUAL VALUE: 0.016ML
- NUMBER OF ACTIVATIONS: 7 DAYS/600 ACTIVATIONS
- UNITS PER BOX: 250

