



PMH-03P-02F REV.01



ONCOFLEX® PRODUCTS

PRODUCTS CATALOG

FL-PMH-014 REV.00 | Data: 20/06/2024

Care with Excellence

PMH - Produtos Médicos Hospitalares, S.A. | +351 263 650 680 | www.pmh.pt | info@pmh.pt



ONCOFLEX® PRODUCTS

EXCELLENCE IN SAFE ADMINISTRATION

ENSURING ASEPTIC PROCEDURES AND MINIMIZING THE RISK OF INFECTIONS IS PMH'S PRIORITY.

THE ONCOFLEX® SYSTEMS ARE CYTOSTATIC ADMINISTRATION/PREPARATION SYSTEMS, DESIGNED TO GUARANTEE MAXIMUM SAFETY AND EFFECTIVENESS IN THE TREATMENT OF PATIENTS.

A VARIANT OF THE ONCOFLEX® SYSTEMS IS MADE UP OF OPAQUE ORANGE TUBES, ESSENTIAL IN PHOTSENSITIVE MEDICATIONS, SERVING AS A BARRIER TO PROTECT PHOTSENSITIVE DRUGS, ENSURING THAT THEY MAINTAIN THEIR INTEGRITY THROUGHOUT THE ADMINISTRATION PROCESS.





ONCOFLEX® SYSTEMS

EXCELLENCE IN SAFE ADMINISTRATION


101.0025 - BLACK OPAQUE SOLUTION ADMINISTRATION SYSTEM

- SPIKE PROTECTION CAP
- SPIKE WITH AERATOR
- ORANGE OPAQUE CHAMBER WITH FILTER 15µm
- WHITE ROLLER CLAMP
- 150CM OF BLACK OPAQUE TUBE
- INJECTION POINT
- FIXED MALE LUER LOCK
- FEMALE WHITE PROTECTIVE CAP
- UNITS PER BOX: 200



101.0026 - ORANGE OPAQUE SOLUTION ADMINISTRATION SYSTEM

- SPIKE PROTECTION CAP
- SPIKE WITH AERATOR
- ORANGE OPAQUE CHAMBER WITH FILTER 15µm
- WHITE ROLLER CLAMP
- 150CM OF ORANGE OPAQUE TUBE
- INJECTION POINT
- FIXED MALE LUER LOCK
- FEMALE WHITE PROTECTIVE CAP
- UNITS PER BOX: 200





ONCOFLEX® SYSTEMS

EXCELLENCE IN SAFE ADMINISTRATION

108.0211 - BLACK OPAQUE SOLUTION ADMINISTRATION SYSTEM WITH LUER

- SPIKE PROTECTION CAP
- SPIKE WITH AERATOR
- ORANGE OPAQUE CHAMBER WITH FILTER 15µm
- WHITE ROLLER CLAMP
- 150CM OF BLACK OPAQUE TUBE
- FIXED MALE LUER LOCK
- FEMALE WHITE PROTECTIVE CAP
- UNIDADES POR CAIXA: 200



108.0212 - ORANGE OPAQUE SOLUTION ADMINISTRATION SYSTEM WITH LUER

- SPIKE PROTECTION CAP
- SPIKE WITH AERATOR
- ORANGE OPAQUE CHAMBER WITH FILTER 15µm
- WHITE ROLLER CLAMP
- 150CM OF ORANGE OPAQUE TUBE
- FIXED MALE LUER LOCK
- FEMALE WHITE PROTECTIVE CAP
- UNIDADES POR CAIXA: 200






ONCOFLEX® SYSTEMS

EXCELLENCE IN SAFE ADMINISTRATION


112.2509 - CYTOSTATIC ADMINISTRATION SYSTEM, WITH 15MM FILTER, MOBILE LUER LOCK AND PURGE FILTER

- SPIKE PROTECTION CAP
- SPIKE WITH AERATOR
- CHAMBER WITH FILTER 15 μ m
- ORANGE ROLLER CLAMP
- FILTER 0.2 μ m
- MOBILE LUER LOCK
- PROTECTIVE CAP WITH PURGE FILTER
- 190CM OF TUBE
- UNITS PER BOX: 100



110.0727 - CYTOSTATIC ADMINISTRATION SYSTEM, WITH 4 PORTS WITH VAR AND MOBILE LUER LOCK AND PURGE FILTER

- SPIKE PROTECTION CAP
- SPIKE WITH AERATOR
- 4 NON-RETURN VALVES (VAR)
- CHAMBER WITH FILTER 15 μ m
- WHITE ROLLER CLAMP
- MOBILE LUER LOCK
- PROTECTIVE CAP WITH PURGE FILTER
- 178CM OF TUBE
- UNITS PER BOX: 40





ONCOFLEX® SYSTEMS

EXCELLENCE IN SAFE ADMINISTRATION

110.1490 - CYTOSTATIC ADMINISTRATION SYSTEM, 2-WAY WITH NEUTROVAL VALVE

- SPIKE PROTECTION CAP
- SPIKE WITH AERATOR
- 2 NEUTROVAL VALVES
- CHAMBER WITH FILTER 15 µm
- ORANGE ROLLER CLAMP
- Y INJECTION POINT WITH NEUTROVAL VALVE
- MOBILE LUER LOCK
- PROTECTIVE CAP WITH PURGE FILTER
- 283CM OF TUBE
- UNITS PER BOX: 40

110.1491 - CYTOSTATIC ADMINISTRATION SYSTEM, 4-WAY WITH NEUTROVAL VALVE




- SPIKE PROTECTION CAP
- SPIKE WITH AERATOR
- 4 NEUTROVAL VALVES
- CHAMBER WITH FILTER 15 µm
- ORANGE ROLLER CLAMP
- Y INJECTION POINT WITH NEUTROVAL VALVE
- MOBILE LUER LOCK
- PROTECTIVE CAP WITH PURGE FILTER
- 216CM OF TUBE
- UNITS PER BOX: 40

112.1334 - OPAQUE, 2-WAY CYTOSTATIC DELIVERY SYSTEM WITH NEUTROVAL VALVE

- SPIKE PROTECTION CAP
- SPIKE WITH AERATOR
- 2 NEUTROVAL VALVES
- CHAMBER WITH FILTER 15 µm
- ORANGE ROLLER CLAMP
- Y INJECTION POINT WITH NEUTROVAL VALVE
- MOBILE LUER LOCK
- PROTECTIVE CAP WITH PURGE FILTER
- 213CM OF ORANGE OPAQUE TUBE
- UNITS PER BOX: 40

112.1338 - OPAQUE, 4-WAY CYTOSTATIC DELIVERY SYSTEM WITH NEUTROVAL VALVE

- SPIKE PROTECTION CAP
- SPIKE WITH AERATOR
- 2 NEUTROVAL VALVES
- CHAMBER WITH FILTER 15 µm
- ORANGE ROLLER CLAMP
- Y INJECTION POINT WITH NEUTROVAL VALVE
- MOBILE LUER LOCK
- PROTECTIVE CAP WITH PURGE FILTER
- 216CM OF ORANGE OPAQUE TUBE
- UNITS PER BOX: 40






ONCOFLEX® EXTENSIONS

EXCELLENCE IN SAFE ADMINISTRATION


313.2348 - CYTOSTATIC PREPARATION LINE WITH SINGLE-CHANNEL SPIKE, BIDIRECTIONAL Y-VALVE

- SPIKE PROTECTION CAP
- SINGLE CHANNEL SPIKE
- Y INJECTION POINT WITH NEUTROVAL VALVE
- MOBILE LUER LOCK
- PROTECTIVE CAP WITH PURGE FILTER
- 38CM OF TUBE
- UNITS PER BOX: 100



313.2346 - CYTOSTATIC PREPARATION LINE WITH SINGLE-CHANNEL SPIKE, BIDIRECTIONAL Y-VALVE

- SPIKE PROTECTION CAP
- SINGLE CHANNEL SPIKE
- Y INJECTION POINT WITH NEUTROVAL VALVE
- TERMINATION WITH NON-RETURN VALVE (VAR)
- PROTECTIVE CAP
- 38CM OF TUBE
- UNITS PER BOX: 100



313.4604 - CYTOSTATIC PREPARATION LINE WITH BICHANNEL SPIKE, Y-SHAPED NEUTROVAL VALVE

- SPIKE PROTECTION CAP
- TWO-CHANNEL SPIKE WITH AERATOR
- Y INJECTION POINT WITH NEUTROVAL VALVE
- FIXED LUER LOCK
- PROTECTIVE CAP WITH PURGE FILTER
- 38CM OF TUBE
- UNITS PER BOX: 100



313.3949 - CYTOSTATIC PREPARATION LINE WITH SINGLE-CHANNEL SPIKE, VAR

- SPIKE PROTECTION CAP
- SINGLE CHANNEL SPIKE
- Y-SHAPED INJECTION POINT WITH NON-RETURN VALVE (VAR)
- FILTER 0.2 µm
- FIXED LUER LOCK
- PROTECTIVE CAP WITH PURGE FILTER
- 38CM OF TUBE
- UNITS PER BOX: 100






ONCOFLEX® EXTENSIONS

EXCELLENCE IN SAFE ADMINISTRATION

313.2405 - CYTOSTATIC PREPARATION LINE WITH TWO-CHANNEL SPIKE, Y-SHAPED NEUTROVAL VALVE AND FILTER

- SPIKE PROTECTION CAP
- TWO-CHANNEL SPIKE WITH AERATOR
- Y INJECTION POINT WITH NEUTROVAL VALVE
- FILTER 0.2 µm
- FIXED LUER LOCK
- PROTECTIVE CAP WITH PURGE FILTER
- 38CM OF TUBE
- UNITS PER BOX: 100



313.1692 - CYTOSTATIC PREPARATION LINE WITH BICHANNEL SPIKE, Y-SHAPED NEUTROVAL VALVE, VAR TERMINATION

- SPIKE PROTECTION CAP
- TWO-CHANNEL SPIKE WITH AERATOR
- Y INJECTION POINT WITH NEUTROVAL VALVE
- FILTER 0.2 µm
- TERMINATION WITH NON-RETURN VALVE (VAR)
- PROTECTIVE CAP
- 38CM OF TUBE
- UNITS PER BOX: 100



313.2430 - OPAQUE CYTOSTATIC PREPARATION LINE WITH Y-SHAPED NEUTROVAL VALVE

- SPIKE PROTECTION CAP
- SINGLE CHANNEL SPIKE
- Y INJECTION POINT WITH NEUTROVAL VALVE
- MOBILE LUER LOCK
- PROTECTIVE CAP WITH PURGE FILTER
- 38CM OF ORANGE OPAQUE TUBE
- UNITS PER BOX: 100



313.4563 - OPAQUE CYTOSTATIC PREPARATION LINE WITH GEN2® Y-VALVE

- SPIKE PROTECTION CAP
- TWO-CHANNEL SPIKE WITH AERATOR
- Y-INJECTION POINT WITH GEN2® VALVE
- MOBILE LUER LOCK
- PROTECTIVE CAP WITH PURGE FILTER
- 38CM OF ORANGE OPAQUE TUBE
- UNITS PER BOX: 100





ONCOFLEX® ACCESSORIES

EXCELLENCE IN SAFE ADMINISTRATION

309.2309 - MINI TWO-CHANNEL SPIKE WITH 0.2 µm FILTER

- SPIKE PROTECTION CAP
- MINI 2-CHANNEL SPIKE
- FILTER 0.2 µm
- LUER LOCK PROTECTION CAP
- UNITS PER BOX: 100



309.0482 - TWO-CHANNEL SPIKE WITH 0.2 µm FILTER

- SPIKE PROTECTION CAP
- BICHANNEL SPIKE
- FILTER 0.2 µm
- LUER LOCK PROTECTION CAP
- UNITS PER BOX: 100



806.2550 - TRANSPARENT SAFETY VALVE

- FEMALE LUER LOCK PROTECTION CAP
- MALE FEMALE LUER LOCK SAFETY VALVE
- PROTECTIVE CAP WITH PURGE FILTER
- UNITS PER BOX: 200



312.2607 - UNIVERSAL ADAPTER FOR SYSTEMS WITH VALVE

- MALE LUER LOCK PROTECTION CAP
- MALE FEMALE LUER LOCK SAFETY VALVE
- 5CM OF TUBE
- PROTECTIVE CAP WITH PURGE FILTER
- UNITS PER BOX: 250





VALVE CHARACTERISTICS

ONCOFLEX® SYSTEM

806.0677 - ANTI-RETURN VALVE (VAR)

- NON-RETURN VALVE.
- LIPID RESISTANT.
- COMPATIBLE WITH CYTOSTATICS AND NMR.
- FLOW ≥ 150 ML/MIN.
- ≥ 90 ML/MIN. (GLUCOSE 40%)
- MAXIMUM PRESSURE RESISTANCE: 6BAR
- OPENING PRESSURE: ≤ 12 MBAR
- UNITS PER BOX: 200



806.2897 - NEUTROVAL® VALVE

- SPLIT-SEPTUM TECHNOLOGY VALVE
- BEFORE EACH USE, THIS MEMBRANE MUST BE DECONTAMINATED WITH 70% ALCOHOL OR CHLORHEXIDINE.
- INJECTION/ASPIRATION WITH DISINFECTION POINT, FOR USE WITHOUT A NEEDLE;
- BARRIER AGAINST MICROORGANISMS IN PROLONGED USE;
- SIMPLE HANDLING AND EASY ACCESS - LUER LOCK AND LUER SLIP;
- COMPATIBLE WITH ALCOHOL, LIPIDS, CHLORHEXIDINE AND CYTOSTATICS;
- COMPATIBLE WITH MRI AND BLOOD.
- PRESSURE TYPE PRESSURE: NEUTRAL
- FLOW RATE: 165ML/MIN
- RESIDUAL VALUE: < 0.04 ML
- MAXIMUM PRESSURE: 27.58 BAR
- NUMBER OF ACTIVATIONS: 7 DAYS / 700 ACTIVATIONS
- UNITS PER BOX: 200



806.4372 - GEN 2® BIDIRECTIONAL VALVE

- PRESSURE TYPE: NEUTRAL
- RESIDUAL VOLUME: 0.016ML
- CONTAINS LATEX: NO
- CONTAINS METAL: NO
- EXTERNAL RESIDUAL FLUID AT DISCONNECTION: MINIMUM
- FIXING SEQUENCE: NOT NECESSARY
- FLOW RATE: 100 ML/MIN
- RETURN PRESSURE: 50PSI
- WEIGHT: 0.6 G
- LENGTH: 26.50 MM
- COMPATIBLE WITH ALCOHOL AND LIPIDS
- UNITS PER BOX: 200

